

TESTS  
FOR  
COLOUR-BLINDNESS

By  
S. ISHIHARA

24 plates

1972



KANEHARA SHUPPAN Co., LTD.

Tokyo • Japan

### Care of the plates

It is important that the book of test plates should be kept closed, except during use, because undue exposure to sunlight causes a fading of the colour of the plates.

Plate	Normal Person	Person with Red-Green Deficiencies		Person with Total Colour Blindness	
		Protan	Deutan	Strong	Mild
1	12	12		12	
2	8	3		×	
3	29	70		×	
4	5	2		×	
5	3	5		×	
6	15	17		×	
7	74	21		×	
8	6	×		×	
9	45	×		×	
10	5	×		×	
11	7	×		×	
12	16	×		×	
13	73	×		×	
14	×	5		×	
15	×	45		×	
		Protan	Deutan		
		Strong	Mild	Strong	Mild
16	26	6	(2) 6	2	2 (6)
17	42	2	(4) 2	4	4 (2)

The mark × shows that the plate cannot be read. Blank space denotes that the reading is indefinite. The numerals in parenthesis show that they can be read but they are comparatively unclear.

































